

Results

Mark Vu offers 14 analysis modes including summary views and deep dives by various skin concerns. Mark Vu can be used as a powerful consultation tool to increase the satisfaction of the client.




QR CODE / PDF FILE SHARE (SKIN ANALYSIS REPORT)

Increase client satisfaction by sharing **personalized skin analysis reports** with QR code and e-mail.




FIXED FOCUS LENS

Mark Vu features a proprietary-developed lens that ensures consistent focus, providing accurate analysis of the skin without any blurring.



HARDWARE ADVANTAGES

It is easy to observe skin problems and compare the before & after with the rotating plate at the bottom. Capturing the left, right, and frontal facial views is convenient with the rotary plate.



Blackout curtain
Completely blocks external light and is easy to use.



One-touch capture button
Simply press the button at the top for easy capturing.

MARK·VU

Professional Facial Skin Analyzer



LIGHT SOURCES

Patent applying No. of light sources 10-2017-0093885
Hidden Wrinkles Application No.10-1981997


Mark-Vu has **full LED light sources**. Observe accurate before & after images and provide customized consultations to your clients.



Normal
Pores
Wrinkles
Skin tone
Dark-Circles



Specular
Hidden Wrinkles
Skin Texture
Radiance




Polarized
Pigmentation
Redness
Brown Spots
Black-Head



UV
Melanin
Sebum
Porphyrin (P.Ance)

CLEAR IMAGES

Easily capture high quality images with our **“Clear View Algorithm.”**



FILTER IN-OUT SYSTEM

Polarizing filter IN & OUT technology enables us to observe the skin with improved color and greater clarity for Normal light and UV light images.



DARK SKIN ANALYSIS

Employing specialized logic, we can detail skin issues in patients with dark skin tones. This is very helpful when capturing images of clients in multi-racial areas **(Europe, South America, the Middle East, etc.)**



SIMULATION MODE

Clients can be shown how **pigmentation, wrinkles, and even hidden wrinkles** may appear without treatment. Consult with our **“Simulation mode”** for effective communication.



RED & BROWN TECHNOLOGY

The **Red & Brown Technology** supports the effective treatment for individuals of diverse ethnic backgrounds.



VIRTUAL LASER TREATMENT

Virtual laser treatment helps to show your clients their expected results of pigmentation treatment and gives them a better understanding of the procedure.



IN-DEPTH ANALYSIS

The **number** and the **ratio** of pores, sebum, and porphyrin can be analyzed with our precise analysis. Pores and pigmentation can be categorized and viewed by three different **size categories** as well.



Size	Number	Average Size
All	4071	0.23 mm ²
Small	2256	0.09 mm ²
Medium	1087	0.20 mm ²
Large	728	0.73 mm ²

Average pore size can be analyzed individually by each size.

3D VIEWER

3D Viewer provides visuals from various dimensions and angles to assess the depth and size of pores, wrinkles, and skin texture.

